

Non-receptor tyrosine-protein kinase TNK1

Catalog No: #AP79777



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Non-receptor tyrosine-protein kinase TNK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD38 negative kinase 1
Accession No.	Uniprot ID:Q13470
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.