

Fibromodulin

Catalog No: #AP79791



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibromodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FM, Collagen-binding 59 kDa protein, Keratan sulfate proteoglycan fibromodulin, KSPG fibromodulin, FM, SLRR2E
Accession No.	Uniprot ID:Q06828
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.