

Fibrinogen-like protein 1

Catalog No: #AP79795



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibrinogen-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mouse fibrinogen-related protein 1,MFIRE-1,Mfire1
Accession No.	Uniprot ID:Q71KU9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.