

Enhancer of filamentation 1

Catalog No: #AP79796



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Enhancer of filamentation 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hEF1,CRK-associated substrate-related protein,CAS-L,Cas scaffolding protein family member 2,Neural precursor cell expressed developmentally down-regulated protein 9,NEDD-9,Renal carcinoma antigen NY-REN-12,p105,CASL,CASS2
Accession No.	Uniprot ID:Q14511
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.