

E3 ubiquitin-protein ligase SIAH2

Catalog No: #AP79800



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase SIAH2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RING-type E3 ubiquitin transferase SIAH2,Seven in absentia homolog 2,Siah-2
Accession No.	Uniprot ID:Q06986
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.