

Kynurenine--oxoglutarate transaminase 1, mitochondrial

Catalog No: #AP79808



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kynurenine--oxoglutarate transaminase 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Kynurenine--oxoglutarate transaminase I, Cysteine-S-conjugate beta-lyase, Glutamine transaminase K, GTK, Glutamine--phenylpyruvate transaminase, Kynurenine aminotransferase 1, Kynurenine aminotransferase I, KATI, Ccb11, Kat
Accession No.	Uniprot ID:Q08415
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.