

Serine/threonine-protein kinase PINK1, mitochondrial

Catalog No: #AP79819



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Serine/threonine-protein kinase PINK1, mitochondrial |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | BRPK,PTEN-induced putative kinase protein 1 |
| Accession No. | Uniprot ID:Q9BXM7 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.