

Vasoactive intestinal polypeptide receptor 1

Catalog No: #AP79825



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vasoactive intestinal polypeptide receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	VIP-R-1,Pituitary adenylate cyclase-activating polypeptide type II receptor,PACAP type II receptor
Accession No.	Uniprot ID:Q28992
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.