

Pancreatic secretory granule membrane major glycoprotein GP2

Catalog No: #AP79827



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pancreatic secretory granule membrane major glycoprotein GP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Pancreatic zymogen granule membrane protein GP-2
Accession No.	Uniprot ID:P25291
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.