

Interferon-induced transmembrane protein 3

Catalog No: #AP79829



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon-induced transmembrane protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Dispanin subfamily A member 2b,DSPA2b
Accession No.	Uniprot ID:P26376
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.