

Interferon-induced transmembrane protein 10

Catalog No: #AP79834



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon-induced transmembrane protein 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dispanin subfamily A member 3,DSPA3
Accession No.	Uniprot ID:Q8BR26
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.