

Probable phosphoglycerate mutase 4

Catalog No: #AP79844



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Probable phosphoglycerate mutase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PGAM3
Accession No.	Uniprot ID:Q8N0Y7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.