

Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma

Catalog No: #AP79855



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Phosphatidylinositol 5-phosphate 4-kinase type II gamma,PI(5)P 4-kinase type II gamma,Pip5k2c |
| Accession No. | Uniprot ID:O88370 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.