

# Protein-L-isoaspartate(D-aspartate) O-methyltransferase

Catalog No: #AP79857



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein-L-isoaspartate(D-aspartate) O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PIMT,L-isoaspartyl protein carboxyl methyltransferase,Protein L-isoaspartyl/D-aspartyl methyltransferase,Protein-beta-aspartate methyltransferase
Accession No.	Uniprot ID:P15246
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.