

Lysine-specific histone demethylase 1A

Catalog No: #AP79861



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lysine-specific histone demethylase 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BRAF35-HDAC complex protein BHC110,Flavin-containing amine oxidase domain-containing protein 2,Aof2,Kiaa0601,Lsd1
Accession No.	Uniprot ID:Q6ZQ88
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.