

Amyloid beta A4 protein

Catalog No: #AP79877



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amyloid beta A4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ABPP,APP,Alzheimer disease amyloid A4 protein homolog,Amyloid precursor protein,Amyloidogenic glycoprotein,AG,Beta-amyloid precursor protein
Accession No.	Uniprot ID:P12023
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.