

6-phosphofructokinase type C

Catalog No: #AP79881



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase type C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	6-phosphofructokinase, platelet type, Phosphofructo-1-kinase isozyme C, PFK-C, Phosphofructokinase 1, Phosphohexokinase, PFKF
Accession No.	Uniprot ID:Q01813
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.