

# Protein arginine N-methyltransferase 1

Catalog No: #AP79884



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein arginine N-methyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Histone-arginine N-methyltransferase PRMT1, Interferon receptor 1-bound protein 4, HMT2, HRMT1L2, IR1B4
Accession No.	Uniprot ID: Q99873
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.