

Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit alpha isoform

Catalog No: #AP79898



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit alpha isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PI3-kinase subunit alpha, Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit alpha, PtdIns-3-kinase subunit p110-alpha, Phosphoinositide-3-kinase catalytic alpha polypeptide, Serine/threonine protein kinase PIK3CA
Accession No.	Uniprot ID:P42336
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.