

Protein Lines homolog 1

Catalog No: #AP79901



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein Lines homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein Lines homolog 2,Wnt-signaling molecule Lines homolog 1,Lins,Lins2,Wins2
Accession No.	Uniprot ID:Q3U1D0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.