

N-lysine methyltransferase KMT5A

Catalog No: #AP79902



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	N-lysine methyltransferase KMT5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H4-K20-HMTase KMT5A, Histone-lysine N-methyltransferase KMT5A, Lysine-specific methylase 5A, PR/SET domain-containing protein 07, PR-Set7, SET domain-containing protein 8, Setd8
Accession No.	Uniprot ID:Q2YDW7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.