

Protein 4.1

Catalog No: #AP79906



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein 4.1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P4.1,4.1R,Band 4.1,Epb4.1,Kiaa4056
Accession No.	Uniprot ID:P48193
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.