

V-set domain-containing T-cell activation inhibitor 1

Catalog No: #AP79924



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	V-set domain-containing T-cell activation inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B7 homolog 4,B7-H4,B7h.5,Immune costimulatory protein B7-H4,Protein B7S1,T-cell costimulatory molecule B7x,B7H4,UNQ659/PRO1291
Accession No.	Uniprot ID:Q7Z7D3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.