

CAAX prenyl protease 1 homolog

Catalog No: #AP79925



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CAAX prenyl protease 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Farnesylated proteins-converting enzyme 1,FACE-1,Prenyl protein-specific endoprotease 1,Zinc metalloproteinase Ste24 homolog,Face1
Accession No.	Uniprot ID:Q80W54
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.