

HERV-K_7p22.1 provirus ancestral Gag polyprotein

Catalog No: #AP79931



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HERV-K_7p22.1 provirus ancestral Gag polyprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EnvK2 protein,Envelope polyprotein,HERV-K(C7) envelope protein,HERV-K(HML-2.HOM) envelope protein,HERV-K108 envelope protein
Accession No.	Uniprot ID:Q69384
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.