

5-hydroxytryptamine receptor 2C

Catalog No: #AP79951



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	5-HT-2C, Serotonin receptor 2C
Accession No.	Uniprot ID:Q60F97
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.