

# 5-hydroxytryptamine receptor 2C

Catalog No: #AP80041



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	5-hydroxytryptamine receptor 2C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5-HT-2C,5-hydroxytryptamine receptor 1C,5-HT-1C,Serotonin receptor 2C,5ht1c,Htr1c
Accession No.	Uniprot ID:P08909
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.