

NT-proBNP

Catalog No: #AP80045



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NT-proBNP
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	N-terminal prohormone of brain natriuretic peptide
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.