

# Aldehyde dehydrogenase family 16 member A1

Catalog No: #AP80046



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Aldehyde dehydrogenase family 16 member A1             |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu                                                     |
| Accession No.      | Uniprot ID:Q8IZ83                                      |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.