

Eukaryotic translation initiation factor 4E-binding protein 1

Catalog No: #AP80053



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eukaryotic translation initiation factor 4E-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4E-BP1, Phosphorylated heat- and acid-stable protein regulated by insulin 1, PHAS-I
Accession No.	Uniprot ID:Q13541
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.