

# Retinal-specific ATP-binding cassette transporter

Catalog No: #AP80067

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Retinal-specific ATP-binding cassette transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATP-binding cassette sub-family A member 4,RIM ABC transporter,RIM protein,Abcr
Accession No.	Uniprot ID:O35600
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.