

4F2 cell-surface antigen heavy chain

Catalog No: #AP80069

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4F2 cell-surface antigen heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	4F2hc, Solute carrier family 3 member 2, CD98
Accession No.	Uniprot ID:Q794F9
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.