

3-hydroxyisobutyrate dehydrogenase, mitochondrial

Catalog No: #AP80072

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-hydroxyisobutyrate dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HIBADH
Accession No.	Uniprot ID:P31937
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.