

3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 1

Catalog No: #AP80088



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type I,3-beta-HSD I,Trophoblast antigen FDO161G,3BH,HSDB3A
Accession No.	Uniprot ID:P14060
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.