

5-hydroxytryptamine receptor 2A

Catalog No: #AP80095



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5-hydroxytryptamine receptor 2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	5-HT-2,Serotonin receptor 2A,5HTR2A
Accession No.	Uniprot ID:Q75Z89
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.