

# Amyloid beta A4 protein

Catalog No: #AP80103



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Amyloid beta A4 protein  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Rt   |
| Other Names        | ABPP,APP,Alzheimer disease amyloid A4 protein homolog,Amyloid precursor protein,Amyloidogenic glycoprotein,AG,Beta-amyloid precursor protein |
| Accession No.      | Uniprot ID:P08592  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.