

Alpha- and gamma-adaptin-binding protein p34

Catalog No: #AP80107



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha- and gamma-adaptin-binding protein p34
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	P34
Accession No.	Uniprot ID:Q9R0Z7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.