

4F2 cell-surface antigen heavy chain

Catalog No: #AP80110



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4F2 cell-surface antigen heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4F2hc,4F2 heavy chain antigen,Lymphocyte activation antigen 4F2 large subunit,Solute carrier family 3 member 2,CD98,MDU1
Accession No.	Uniprot ID:P08195
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.