

Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial

Catalog No: #AP80136

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	KAT/AadAT, 2-aminoadipate aminotransferase, 2-aminoadipate transaminase, Alpha-aminoadipate aminotransferase, AadAT, Kynurenine aminotransferase II, Kynurenine--oxoglutarate aminotransferase II, Kynurenine--oxoglutarate transaminase 2, Kynurenine--oxoglutarate transaminase II
Accession No.	Uniprot ID: Q5E9N4
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.