

Aspartate aminotransferase, mitochondrial

Catalog No: #AP80150

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aspartate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mAspAT,Fatty acid-binding protein,FABP-1,Glutamate oxaloacetate transaminase 2,Kynurenine aminotransferase 4,Kynurenine aminotransferase IV,Kynurenine--oxoglutarate transaminase 4,Kynurenine--oxoglutarate transaminase IV,Plasma membrane-associated fatty acid-binding protein,FABPpm,Transaminase A,Got-2
Accession No.	Uniprot ID:P05202
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.