

# Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial

Catalog No: #AP80162

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	KAT/AadAT,2-aminoadipate aminotransferase,2-aminoadipate transaminase,Alpha-aminoadipate aminotransferase,AadAT,Kynurenine aminotransferase II,Kynurenine--oxoglutarate aminotransferase II,Kynurenine--oxoglutarate transaminase 2,Kynurenine--oxoglutarate transaminase II,Kat2
Accession No.	Uniprot ID:Q9WVM8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.