

Long-chain-fatty-acid--CoA ligase 1

Catalog No: #AP80185



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Long-chain-fatty-acid--CoA ligase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acyl-CoA synthetase 1,ACS1,Long-chain acyl-CoA synthetase 1,LACS 1,Long-chain acyl-CoA synthetase 2,LACS 2,Long-chain fatty acid-CoA ligase 2,Palmitoyl-CoA ligase 1,Palmitoyl-CoA ligase 2,FACL1,FACL2,LACS,LACS1,LACS2
Accession No.	Uniprot ID:P33121
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.