

5'-AMP-activated protein kinase catalytic subunit alpha-1

Catalog No: #AP80233



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5'-AMP-activated protein kinase catalytic subunit alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	AMPK subunit alpha-1,AMPK 63 kDa subunit,Acetyl-CoA carboxylase kinase,ACACA kinase,Hydroxymethylglutaryl-CoA reductase kinase,HMGCR kinase,Tau-protein kinase PRKAA1,AMPK1
Accession No.	Uniprot ID:Q09136
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.