

## A-kinase anchor protein 1, mitochondrial

Catalog No: #AP80263

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A-kinase anchor protein 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	A-kinase anchor protein 121 kDa,AKAP 121,Dual specificity A-kinase-anchoring protein 1,D-AKAP-1,Protein kinase A-anchoring protein 1,PRKA1,Spermatid A-kinase anchor protein 84,S-AKAP84,Akap
Accession No.	Uniprot ID:O88884
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.