

# Serum albumin

Catalog No: #AP80328



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serum albumin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GIG20,GIG42,PRO0903,PRO1708,PRO2044,PRO2619,PRO2675,UNQ696/PRO1341
Accession No.	Uniprot ID:P02768
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.