

Serine/threonine-protein kinase receptor R3

Catalog No: #AP80336



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase receptor R3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SKR3,Activin receptor-like kinase 1,ALK-1,TGF-B superfamily receptor type I,TSR-I,Acvr1k1,Alk-1
Accession No.	Uniprot ID:Q61288
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.