

# Aryl hydrocarbon receptor

Catalog No: #AP80367



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Aryl hydrocarbon receptor  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Ah receptor, Class E basic helix-loop-helix protein 76, bHLHe76, BHLHE76 |
| Accession No.      | Uniprot ID:P35869  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                   |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.