

Activity-dependent neuroprotector homeobox protein

Catalog No: #AP80381



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Activity-dependent neuroprotector homeobox protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activity-dependent neuroprotective protein
Accession No.	Uniprot ID:Q9JKL8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.