

# Protein argonaute-2

Catalog No: #AP80384



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein argonaute-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Argonaute2,Argonaute RISC catalytic component 2,Eukaryotic translation initiation factor 2C 2,eIF-2C 2,Piwi/argonaute family protein melf2C2,Protein slicer,Eif2c2,Kiaa4215
Accession No.	Uniprot ID:Q8CJG0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.