

Type-1A angiotensin II receptor

Catalog No: #AP80385



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Type-1A angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin II type-1A receptor,AT1A,Agtr1
Accession No.	Uniprot ID:P29754
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.